

HOSPITAL



# LIFETHERM

Warming beds for preterm infants, neonates and babies

## LIFETHERM 2000

### Warming bed

The warming bed is available without cover, with a closed cover or with a fold up cover that can be opened on both sides, and can be upgraded with extensive accessories. It provides even heat distribution due to the heated gel mattress, which is an ideal area for the neonate. Lifetherm 2000 is the perfect solution for the intensive care of the small patients and is characterized by easy usability.

- Effective heat transfer due to panel heater
- Precision electronic adjustment of set temperature
- Display of current temperature
- Infinite electronic height adjustment with footswitches on both sides\*
- Gentle electric bed tilting by pressing a button
- Front and side panels are prepared to turn down
- Covers available that can be opened on either side, as well as one-piece covers\*
- Gel mattress



## lifetherm 2000



\* optional accessory

## LIFETHERM 2002 / 2003

### Resuscitation unit

A family of resuscitation unit with a lot of different options for optimal adjustment to suit the application as well as the needs of the baby and operating staff. Lifetherm 2002 and Lifetherm 2003 create an excellent heatarea for the small patient.

- Swivelling, fast-response radiant heater with integrated LED examination lamp
- Conductive heat by a heated gel mattress<sup>1</sup>
- Quiet, electric height adjustment with footswitch control on both sides of the bed\*
- Gentle electric bed tilting by pressing the corresponding button
- Fold-down side panels available in various heights
- Skin temperature control and display (accuracy  $\pm 0.1^{\circ}\text{C}$ )\*\*
- Covers available that can be opened on either side, as well as one-piece covers\*
- Trays which allow the attachment of control panels and ventilators
- Standard rails for attaching accessories such as catheter holders, trays or suction devices\*



## lifetherm 2002 / 2003

### Benefits of the Lifetherm series at a glance:

- Generated contact warmth by gel mattress
- Large-area radiant warmth due to ceramic lamp\*
- Ergonomic device with easy-clean design
- Manual skin temperature control\*
- Lockable casters
- Electrical height adjustment with control from both sides of the bed\*
- Gentle electric bed adjustment for head-up and head-down positions
- Front and side panels are prepared to turn down
- Wide range of accessories

A main priority of gentle care is to reduce the level of stress for the small patient, its relatives, and of course nursing staff. In conjunction with the experience of Löwenstein Medical in the field of warming therapy this realization has led to a user-friendly, development-oriented family of equipment. An extensive range of accessories allows it to adapt to particular specific requirements.

\* optional accessory

\*\* only Lifetherm 2002

<sup>1</sup> Gel mattress only in conjunction with the heated bed surface



Löwenstein Medical  
Arzbacher Strasse 80  
D-56130 Bad Ems

Tel.: +49 26 03/96 00-0  
Fax: +49 26 03/96 00-50  
Internet: hul.de

*The images of our products contain optional accessories which can be ordered in addition to the Product.  
Please contact for more information.*

to the product category:

